

NIR Spectrometers 1000 – 1650 nm

(low cost, high sensitivity, USB)

Patent pending



DATASHEET

[Return to the Webpage](#)



Features

- High Sensitivity
- Low Cost
- USB/GUI
- Deep-Cooling Option
- MEMS Chopper option

Applications

- Laboratory use
- Testing
- Instrumentation

The NIRS Series Spectrometer, based on TE-cooled InGaAs detector coupled with a grating, is optimized for NIR spectroscopic measurements requiring an exceptional signal-to-noise ratio and high dynamic range across the 0.9 - 2.5 μm spectral range. This system leverages innovative, patent-pending scanning technology, offering significant advantages: 1) Unmatched low cost; 2) Industry-leading sensitivity with deep cooling to -40°C ; 3) Extended spectral coverage beyond traditional spectrometers; 4) Low power consumption; 5) Integrated MEMS chopper; 6) High- resolution performance. Additional features include photon integration for low-noise detection and connectivity via USB or RS232 with an intuitive GUI. The NIRS Series spectrometers deliver high performance with ultra-low noise levels, making them suitable for a range of demanding applications. The detectors' excellent sensitivity supports broad-band applications, such as analyzing the optical properties of solids, liquids, and gases in the NIR range, chemical component analysis, moisture detection, and narrow-bandwidth tasks like NIR laser characterization. The NIRS series comes standard with a USB interface, and software support includes SDK examples, DLLs for custom application development, and Windows-based spectral acquisition and analysis tools.

Specifications

Parameter		Min	Typical	Max	Unit
Wavelength Range		1000		1650	nm
Resolution		0.2	0.4	0.6	nm
Wavelength Accuracy		0.05	0.08	0.1	nm
PDL		-	0.15	0.35	dB
Sensor Size			0.2		mm
Noise Floor ^[1]		-110		-10	dBm
Scan Time		1		200	s
Input Optical Power	Standard version	-		0.3	W
	High power version			5	W
f-number			1.8		
Electronic Interface				Mini USB	
Operating Temperature		0	20	60	$^{\circ}\text{C}$
Storage Temperature		-10	-	70	$^{\circ}\text{C}$
Weight			10		lbs

Notes:

[1]. The lowest level requires -40°C cooling, the high level is room temperature. These are also related to the integration time setting. Low spectral resolution increase sensitivity.

[2]. An integrated shutter is available to calibrate the dark readout

Note: The specifications provided are for general applications with a cost-effective approach. If you need to narrow or expand the tolerance, coverage, limit, or qualifications, please [click this link](#):

Legal notices: All product information is believed to be accurate and is subject to change without notice. Information contained herein shall legally bind Agiltron only if it is specifically incorporated into the terms and conditions of a sales agreement. Some specific combinations of options may not be available. The user assumes all risks and liability whatsoever in connection with the use of a product or its application.

Rev 05/08/25

© Photonwares Corporation

P +1 781-935-1200

E sales@photonwares.com

W www.agiltron.com

Information contained herein is deemed to be reliable and accurate as of the issue date. Photonwares reserves the right to change the design or specifications at any time without notice. Agiltron is a registered trademark of Photonwares Corporation in the U.S. and other countries.

NIR Spectrometers 1000 – 1650 nm


(low cost, high sensitivity, USB)

Patent pending



DATASHEET

Dimensions (mm)



*Product dimensions may change without notice. This is sometimes required for non-standard specifications.

Electrical/Computer Connection

Module comes with a 12V DC power wall pluggable power supply

Turn-key unit has a 110-220 ACV input and a USB input at the back and optical input at the front.

USB 2.0, 480 Mbps; RS-232 (3-wire)

Ordering Information

	1	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Prefix	Type	Wavelength	Input Optical Power	Cooling *	Resolution **	Shutter	Chopper	Connector
NIRS-		1-1.6 μ m = 1 Special = 0	Standard = 1 High Power = 2	Non = 1 -5°C = 2 -20°C = 3 -40°C = 5	1nm = 1 0.5nm = 2 5nm = 3 10nm = 4	Non = 1 Yes = 2	Non = 1 Yes = 2	SMA905 = 1 FC/PC = 2 SC/PC = 4 ST/PC = 6 Special = 0

* Non cooling is low cost for strong light measurements. At -5°C: Noise is reduced by about 4 \times , improving performance in low-light applications. At -40°C: Noise is reduced by about 16 \times , enabling high-sensitivity measurements, such as weak signal detection in spectroscopy or astronomy.

** Low resolution high sensitivity.